KULLANIM AMACI

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INTENDED USE

(English)

Diachek 3D is a hematology control tested to evaluate the accuracy and precision of hematology devices providing white blood cell differential.

SUMMARY AND PRINCIPLES

Laboratories require assayed material for quality control of automated, semi-automatic and manual procedures measuring whole blood parameters. Daily use of this whole blood control provides quality control data to confirm the precision and accuracy of device operation.

The use of stabilized cell preparations to control hematology analyzers is an established procedure. The Diachek 3 DIFF, when handled like a patient sample and tested on a properly calibrated and working instrument, will provide values within the expected range specified on the test sheet.

CONTROL MATERIAL

This whole blood control material may contain any or all of the following: A platelet component consisting of stabilized human red blood cells, a white blood cell component composed of human and/or non-human analogues, and a human or non-human analogue in a preservative.

PRECAUTIONS

1. For IVD use only.
2. Although the serology tests prescribed by WHO were performed and negative samples were used, it should not be forgotten that it is a blood material and may contain infectious virus.
3. It should be treated as medical waste and disposed of.

STORAGE AND STABILITY

It should be stored at 2-8 °C. Stable from unopened to expiration date, 20 days after capping.

PRODUCT DEFAULT

It may be an indication that it is outside the expected values. However, excessive darkening in the color of the product or formation of clots indicates that the product has deteriorated.

USAGE

1. Remove the control vial at 2-8 °C and leave at room temperature for 15 minutes.
2. You can mix it in two different ways.
 - a. Hold the bottle vertically and rotate each bottle between your palms for 15-20 seconds.

- b. Holding the vial at the ends between thumb and finger, continue mixing by rapidly inverting the vial 20

times from end to end using a very quick whirling motion of the wrist.

- c. Analyze immediately after mixing. Device analysis can be performed by inverting the vial 5 times.
3. See instrument manual for system used to analyze control materials.
 4. After use, the bottle should be closed tightly and immediately put in the refrigerator.

LIMITATIONS

Manual differential analysis of white blood cells cannot be performed with DiaCheck 3 DIFF. White blood cell components simulate white blood cells in size, not morphology.

EXPECTED RESULTS

The average assay values provided for each parameter are derived from repeated analysis on calibrated instruments as well as manual reference methods. MCV and hematocrit values are not corrected for entrapped plasma. Assay values are obtained using reagents recommended by instrument manufacturers and will be used for instrument control; these are not absolute tests for calibration.

When a new lot of control is purchased, it is good laboratory practice for a separate laboratory to set its own mean and limits for each parameter. However, the control means created by the laboratory should be within the expected range specified for the control. The expected ranges listed represent estimates of variation due to different laboratories, instrument calibration, maintenance and operator technique..

For devices not listed on the test page, test values and expected ranges must be determined by the user. In accordance with CLSI H26-A2¹, it is recommended to run each control level twice a day for 3-5 days to create separate laboratory instruments for each measurand.

REFERENCES

1. Clinical and Laboratory Standards Institute, H26-A2, Validation, verification, and quality assurance of automated hematology analyzers; Approved Standard – Second Edition.

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Medios Medikal Bilişim Danışmanlık Tic. Ltd. Şti.